

APLW 036

Lid

height 1 49/64" (45 mm), for fixing delicate wash items.



EAN: 4002516115984 / Material number: 11104960 / Old Material Number: 69303600D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Other inserts	•
Cover	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in inches	9 13/16 (248)
External dimensions, net width in inches	9 13/16 (248)
External dimensions, net depth in inches	1 13/16 (45)
External dimensions, gross height in inches	11 (278)
External dimensions, gross width in inches	2 1/4 (57)
External dimensions, gross depth in inches	11 1/4 (285)
Net weight in lbs	1 (0.39)
Gross weight in lbs	1 (0.54)