

APLW 039
 Mesh Basket
 for optimum loading of half petri dishes.



- Suitable for 26 half petri dishes
- Distance between the holders = 15/16" (24 mm.)

EAN: 4002516116011 / Material number: 11104990 / Old Material Number: 69303900D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for half petri dishes	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	3 5/16 (83)
External dimensions, net width in inches	9 3/8 (237)
External dimensions, net depth in inches	16 15/16 (430)
External dimensions, gross height in inches	10 1/16 (255)
External dimensions, gross width in inches	16 3/16 (410)
External dimensions, gross depth in inches	20 1/8 (510)
Net weight in lbs	3 (1.17)
Gross weight in lbs	6 (2.47)