


APLW 040  
 Mesh Basket  
 for optimum loading of test tubes, height 3 15/16" (100 mm.)

EAN: 4002516116028 / Material number: 11105000 / Old Material Number: 69304000D



- For use with APLW 043 lid

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
Maximum vertical clearance in mm	4 9/16 (115)
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	3 15/16 (100)
External dimensions, net width in inches	8 11/16 (220)
External dimensions, net depth in inches	8 11/16 (220)
External dimensions, gross height in inches	13 (330)
External dimensions, gross width in inches	12 7/16 (315)
External dimensions, gross depth in inches	12 7/16 (315)
Net weight in lbs	2 (0.47)
Gross weight in lbs	4 (1.58)