


APLW 041  
 Mesh Basket  
 for optimum loading of test tubes, height 5 1/8" (130 mm.)

EAN: 4002516116035 / Material number: 11105010 / Old Material Number: 69304100D



- For use with APLW 043 lid

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
Maximum vertical clearance in mm	4 5/8 (115)
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	130
External dimensions, net width in inches	220
External dimensions, net depth in inches	220
External dimensions, gross height in inches	13 (330)
External dimensions, gross width in inches	11 7/8 (300)
External dimensions, gross depth in inches	12 1/2 (315)
Net weight in lbs	2 (0.55)
Gross weight in lbs	3 (1.35)