


APLW 043
Lid
for APLW 040-042 inserts



EAN: 4002516116059 / Material number: 11105030 / Old Material Number: 69304300D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Lid	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	0 5/8 (15)
External dimensions, net width in inches	8 1/8 (205)
External dimensions, net depth in inches	8 1/8 (205)
External dimensions, gross height in inches	8 1/2 (215)
External dimensions, gross width in inches	6 3/4 (170)
External dimensions, gross depth in inches	10 11/16 (270)
Net weight in lbs	1 (0.18)
Gross weight in lbs	2 (0.55)