


APLW 044
 Partition mesh
 for dividing the base area in APLW 040–042.



EAN: 4002516116066 / Material number: 11105040 / Old Material Number: 69304400D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Insert/module for test tubes	•
Holder	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	3 (75)
External dimensions, net width in inches	8 1/8 (205)
External dimensions, net depth in inches	7 7/8 (200)
External dimensions, gross height in inches	8 1/2 (215)
External dimensions, gross width in inches	6 3/4 (170)
External dimensions, gross depth in inches	10 11/16 (270)
Net weight in lbs	1 (0.36)
Gross weight in lbs	2 (0.60)