

APLW 047

Injector nozzle

with plastic support, conical, Ø 3/32" (2.5 mm), length 2" (50 mm.)

EAN: 4002516117261 / Material number: 11105070 / Old Material Number: 69304700D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 5/32" (4 mm.)

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	1 5/16 (32)
External dimensions, net width in inches	1 5/16 (32)
External dimensions, net depth in inches	2 7/16 (61)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.02)