

APLW 050

Injector nozzle

with plastic support, star-shaped, Ø 5/32" (4 mm), length 2 61/64" (75 mm.)

EAN: 4002516116110 / Material number: 11105100 / Old Material Number: 69305000D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 13/64" (5 mm.)

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	1 5/16 (33)
External dimensions, net width in inches	1 1/2 (37)
External dimensions, net depth in inches	3 1/16 (77)
External dimensions, gross height in inches	1 5/8 (40)
External dimensions, gross width in inches	2 9/16 (65)
External dimensions, gross depth in inches	4 11/16 (118)
Net weight in lbs	1 (0.01)
Gross weight in lbs	1 (0.03)