


APLW 051

Injector nozzle

with plastic support, star-shaped with locking retainer, Ø 4, length 110 mm.

EAN: 4002516116127 / Material number: 11105110 / Old Material Number: 69305100D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 13/64" (5 mm.)

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	2 3/16 (54) in (mm)
External dimensions, net width in inches	2 3/16 (54) in (mm)
External dimensions, net depth in inches	4 9/16 (115) in (mm)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.01)
Gross weight in lbs	1 (0.02)