

APLW 052

Injector nozzle

with plastic support, star-shaped with locking retainer, Ø 4, length 175 mm.

EAN: 4002516117278 / Material number: 11105120 / Old Material Number: 69305200D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 13/64" (5 mm.)

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	2 3/16 (54)
External dimensions, net width in inches	2 3/16 (54)
External dimensions, net depth in inches	7 1/16 (178)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 15/16 (100)
External dimensions, gross depth in inches	5 15/16 (150)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.02)