

APLW 055

Injector nozzle

with plastic support, star-shaped with locking retainer, Ø 6, length 175 mm.

EAN: 4002516116158 / Material number: 11105150 / Old Material Number: 69305500D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 25/64" (10 mm.)

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	3 (75)
External dimensions, net width in inches	3 (75)
External dimensions, net depth in inches	7 3/16 (182)
External dimensions, gross height in inches	2 3/16 (55)
External dimensions, gross width in inches	7 1/8 (180)
External dimensions, gross depth in inches	9 7/8 (250)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.03)