

APLW 059

Injector nozzle

with plastic support & bottle support with locking retainer, Ø 6, length 225 mm.

EAN: 4002516116257 / Material number: 11105190 / Old Material Number: 69305900D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 25/64" (10 mm.)

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	3 1/16 (77)
External dimensions, net width in inches	3 1/16 (77)
External dimensions, net depth in inches	10 3/16 (258)
External dimensions, gross height in inches	2 3/16 (55)
External dimensions, gross width in inches	3 15/16 (100)
External dimensions, gross depth in inches	11 13/16 (300)
Net weight in lbs	0 (0.06)
Gross weight in lbs	0 (0.12)