

APLW 062

Injector nozzle

with plastic support and bottle neck holder, star-shaped, Ø 6, length 135 mm.

EAN: 4002516116288 / Material number: 11105220 / Old Material Number: 69306200D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 25/64" (10 mm.)
- With locking retainer

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	76
External dimensions, net width in inches	76
External dimensions, net depth in inches	145
External dimensions, gross height in inches	3 1/4 (80)
External dimensions, gross width in inches	3 3/4 (95)
External dimensions, gross depth in inches	4 3/8 (110)
Net weight in lbs	1 (0.03)
Gross weight in lbs	1 (0.07)