

APLW 064

Injector nozzle

with plastic support and bottle neck holder, star-shaped, Ø 6, length 225 mm.

EAN: 4002516116301 / Material number: 11105240 / Old Material Number: 69306400D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 25/64" (10 mm.)
- With locking retainer

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	3 (76)
External dimensions, net width in inches	3 (76)
External dimensions, net depth in inches	9 1/8 (231)
External dimensions, gross height in inches	2 3/16 (55)
External dimensions, gross width in inches	3 15/16 (100)
External dimensions, gross depth in inches	11 13/16 (300)
Net weight in lbs	1 (0.04)
Gross weight in lbs	1 (0.10)