

APLW 065

Injector nozzle

with plastic support and bottle neck holder, star-shaped, Ø 6, length 275 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 25/64" (10 mm.)
- With locking retainer

EAN: 4002516116318 / Material number: 11105250 / Old Material Number: 69306500D

|                                                |                          |
|------------------------------------------------|--------------------------|
| <b>Product category</b>                        |                          |
| Large-capacity washer-disinfectors, laboratory | •                        |
| <b>Machine affiliation</b>                     |                          |
| PLW 6011                                       | •                        |
| PLW 6111                                       | •                        |
| <b>Product category</b>                        |                          |
| Injector nozzle                                | •                        |
| <b>Product characteristics</b>                 |                          |
| Material                                       | Stainless steel, Plastic |
| Color                                          | Blue, Stainless Steel    |
| <b>Application area</b>                        |                          |
| Reprocessing of laboratory glassware           | •                        |
| <b>Dimensions and weight</b>                   |                          |
| External dimensions, net height in inches      | 3 (75)                   |
| External dimensions, net width in inches       | 3 (75)                   |
| External dimensions, net depth in inches       | 11 5/16 (287)            |
| External dimensions, gross height in inches    | 2 3/16 (55)              |
| External dimensions, gross width in inches     | 3 15/16 (100)            |
| External dimensions, gross depth in inches     | 11 13/16 (300)           |
| Net weight in lbs                              | 1 (0.04)                 |
| Gross weight in lbs                            | 1 (0.11)                 |