

APLW 068

Injector nozzle

with holder, Ø 15/64" (6 mm), length 10 53/64" (275 mm.)



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 25/64" (10 mm.)

EAN: 4002516116349 / Material number: 11105280 / Old Material Number: 69306800D

| Product category | |
|--|--------------------------------------|
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector nozzle | • |
| Product characteristics | |
| Material | Wire, rilsan-coated, Stainless steel |
| Color | Blue, Stainless Steel, White |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in inches | 70 |
| External dimensions, net width in inches | 85 |
| External dimensions, net depth in inches | 287 |
| External dimensions, gross height in inches | 2 1/4 (55) |
| External dimensions, gross width in inches | 4 (100) |
| External dimensions, gross depth in inches | 11 7/8 (300) |
| Net weight in lbs | 1 (0.09) |
| Gross weight in lbs | 1 (0.15) |