

APLW 068
 Injector nozzle
 with holder, Ø 15/64" (6 mm), length 10 53/64" (275 mm.)

EAN: 4002516116349 / Material number: 11105280 / Old Material Number: 69306800D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Including nozzle cap, Ø 25/64" (10 mm.)

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Wire, rilsan-coated, Stainless steel
Color	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	2 13/16 (70)
External dimensions, net width in inches	3 3/8 (85)
External dimensions, net depth in inches	11 5/16 (287)
External dimensions, gross height in inches	2 3/16 (55)
External dimensions, gross width in inches	3 15/16 (100)
External dimensions, gross depth in inches	11 13/16 (300)
Net weight in lbs	1 (0.09)
Gross weight in lbs	1 (0.15)