

APLW 069

Injector nozzle

with stainless steel support, Ø 5/16" (8 mm), length 10 3/64" (255 mm.)

EAN: 4002516116356 / Material number: 11105290 / Old Material Number: 69306900D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Grey
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	3 7/16 (87)
External dimensions, net width in inches	3 7/16 (87)
External dimensions, net depth in inches	10 3/8 (263)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	9 7/8 (250)
External dimensions, gross depth in inches	13 13/16 (350)
Net weight in lbs	1 (0.08)
Gross weight in lbs	1 (0.09)