

APLW 070

Injector nozzle

with stainless steel support, Ø 5/16" (8 mm), length 12 19/32" (320 mm.)



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

EAN: 4002516116363 / Material number: 11105300 / Old Material Number: 69307000D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Grey
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	87
External dimensions, net width in inches	87
External dimensions, net depth in inches	330
External dimensions, gross height in inches	4 3/8 (108)
External dimensions, gross width in inches	4 5/8 (117)
External dimensions, gross depth in inches	12 1/8 (305)
Net weight in lbs	1 (0.10)
Gross weight in lbs	1 (0.19)