

APLW 070

Injector nozzle

with stainless steel support, Ø 5/16" (8 mm), length 12 19/32" (320 mm.)



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

EAN: 4002516116363 / Material number: 11105300 / Old Material Number: 69307000D

| | |
|--|--------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector nozzle | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Color | Stainless Steel, Grey |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in inches | 3 7/16 (87) |
| External dimensions, net width in inches | 3 7/16 (87) |
| External dimensions, net depth in inches | 13 (330) |
| External dimensions, gross height in inches | 4 5/16 (108) |
| External dimensions, gross width in inches | 4 5/8 (117) |
| External dimensions, gross depth in inches | 12 1/16 (305) |
| Net weight in lbs | 1 (0.10) |
| Gross weight in lbs | 1 (0.19) |