

APLW 071

Injector nozzle

w/ holder for meas. cylinders up to Ø 105 mm & inj. nozzles Ø 8, length 320 mm.

EAN: 4002516116370 / Material number: 11105310 / Old Material Number: 69307100D



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring cylinders, flasks, and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Grey
<b>Dimensions and weight</b>	
External dimensions, net height in inches	4 3/4 (120)
External dimensions, net width in inches	4 3/4 (120)
External dimensions, net depth in inches	13 (330)
External dimensions, gross height in inches	4 15/16 (125)
External dimensions, gross width in inches	6 1/2 (165)
External dimensions, gross depth in inches	13 13/16 (350)
Net weight in lbs	1 (0.24)
Gross weight in lbs	1 (0.39)