

APLW 086
Bottle neck holder
Ø 1 49/64" (45 mm), for injector nozzles Ø 15/64" (6 mm.)



EAN: 4002516116646 / Material number: 11105460 / Old Material Number: 69308600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other accessories	•
Product characteristics	
Material	Plastic
Color	White
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	1 5/16 (33)
External dimensions, net width in inches	1 5/16 (33)
External dimensions, net depth in inches	2 13/16 (70)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	4 3/4 (120)
External dimensions, gross depth in inches	6 3/4 (170)
Net weight in lbs	1 (0.05)
Gross weight in lbs	1 (0.05)