

APLW 088

Bracket

for measuring cylinders up to Ø 105 mm and injector nozzles Ø 8.



EAN: 4002516116660 / Material number: 11105480 / Old Material Number: 69308800D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Holder	•
Other accessories	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	12
External dimensions, net width in inches	120
External dimensions, net depth in inches	300
External dimensions, gross height in inches	4 3/8 (108)
External dimensions, gross width in inches	4 5/8 (117)
External dimensions, gross depth in inches	12 1/8 (305)
Net weight in lbs	1 (0.14)
Gross weight in lbs	1 (0.23)