

APLW 090

Holder

for pipettes with max. 35/64" (14 mm) and min. 5/32" (4 mm) diameter.



EAN: 4002516116684 / Material number: 11105500 / Old Material Number: 69309000D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Holder	•
Other accessories	•
<b>Product characteristics</b>	
Material	Plastic
Color	Grey, Black
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	0 13/16 (20)
External dimensions, net width in inches	0 13/16 (20)
External dimensions, net depth in inches	0 13/16 (20)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	4 3/4 (120)
External dimensions, gross depth in inches	6 3/4 (170)
Net weight in lbs	1 (0.00)
Gross weight in lbs	1 (0.00)