

APLW 091
Holder
for injector nozzles Ø 3/32" (2.5 mm.)



EAN: 4002516116691 / Material number: 11105510 / Old Material Number: 69309100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Holder	•
Other accessories	•
Product characteristics	
Material	Plastic
Color	Grey
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	0 13/16 (20)
External dimensions, net width in inches	0 13/16 (20)
External dimensions, net depth in inches	0 11/16 (17)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.00)
Gross weight in lbs	1 (0.00)