

APLW 093  
 Irrigation sleeves  
 with cap for pipettes, max. Ø 7/16" (11 mm.)



EAN: 4002516116721 / Material number: 11105530 / Old Material Number: 69309300D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Holder	•
Other accessories	•
<b>Product characteristics</b>	
Material	Stainless steel, Silicone
Color	Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	1 (25)
External dimensions, net width in inches	1 (25)
External dimensions, net depth in inches	4 11/16 (118)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.05)
Gross weight in lbs	1 (0.05)