

A 841/3

Injector nozzle

With plastic supports, nozzle diameter 6, length 210 mm, 20 items



- Suitable for mobile injector units & baskets, w/ screw connection
- Part of the SimpleLoad system
- Packaging unit with 20 nozzles

EAN: 4002516188834 / Material number: 11113150 / Old Material Number: 69684103D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	8 5/16 (210)
External dimensions, net width in inches	3 5/8 (92)
External dimensions, net depth in inches	3 7/16 (86)
External dimensions, gross height in inches	5 9/16 (140)
External dimensions, gross width in inches	7 1/2 (190)
External dimensions, gross depth in inches	11 5/8 (295)
Net weight in lbs	2 (0.73)
Gross weight in lbs	2 (0.89)