

A 840/2

Injector nozzle

With plastic supports, nozzle diameter 6, length 130 mm, 10 items.



- Suitable for mobile injector units & baskets, w/ screw connection
- Part of the SimpleLoad system
- Packaging unit with 10 nozzles

EAN: 4002516188841 / Material number: 11113160 / Old Material Number: 69684002D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	5 1/8 (130)
External dimensions, net width in inches	3 5/8 (92)
External dimensions, net depth in inches	3 7/16 (86)
External dimensions, gross height in inches	5 9/16 (140)
External dimensions, gross width in inches	5 15/16 (150)
External dimensions, gross depth in inches	10 1/4 (260)
Net weight in lbs	1 (0.27)
Gross weight in lbs	1 (0.40)