

A 843/1

Injector nozzle

With plastic supports, nozzle diameter 4, length 185 mm, 5 items



- Suitable for mobile injector units & baskets, w/ screw connection
- Part of the SimpleLoad system
- Packaging unit with 5 nozzles

EAN: 4002516188889 / Material number: 11113720 / Old Material Number: 69684301D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	7 5/16 (185)
External dimensions, net width in inches	3 (76)
External dimensions, net depth in inches	2 13/16 (71)
External dimensions, gross height in inches	4 1/4 (107)
External dimensions, gross width in inches	5 1/8 (130)
External dimensions, gross depth in inches	7 11/16 (195)
Net weight in lbs	1 (0.11)
Gross weight in lbs	1 (0.22)