

A 845/1

Injector nozzle

With plastic supports, nozzle diameter 2.5, length 125 mm, 5 items.



- Suitable for mobile injector units & baskets, w/ screw connection
- Part of the SimpleLoad system

EAN: 4002516188940 / Material number: 11116120 / Old Material Number: 69684501D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 7111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	4 15/16 (125)
External dimensions, net width in inches	2 3/8 (60)
External dimensions, net depth in inches	2 3/16 (55)
External dimensions, gross height in inches	3 3/8 (85)
External dimensions, gross width in inches	4 3/4 (120)
External dimensions, gross depth in inches	5 15/16 (150)
Net weight in lbs	1 (0.06)
Gross weight in lbs	1 (0.13)