

A 300/3
Module

For optimum loading of laboratory glassware – can be equipped with 8 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Equipped with 4 A 840 nozzles and 4 A 841 nozzles
- Part of the SimpleLoad system
- Nozzle spacing 4 5/8" across the width, 4 5/8" across the depth

EAN: 4002516188964 / Material number: 11116900 / Old Material Number: 69630003D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Number of compartments	8
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	8
Dimensions and weight	
External dimensions, net height in inches	9 (228) in (mm)
External dimensions, net width in inches	8 1/4 (208) in (mm)
External dimensions, net depth in inches	18 7/8 (479) in (mm)
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	9 1/2 (240)
External dimensions, gross depth in inches	20 5/16 (515)
Net weight in lbs	4 (1.47)
Gross weight in lbs	5 (1.91)
Accessory pack	
A 840, injector nozzle Ø 1/4" (6 mm), length 5 1/8" (130 mm) [quantity]	4
A 841, injector nozzle Ø 1/4" (6 mm), length 8 1/4" (210 mm) [quantity]	4