


A 301/5

Module

For optimum loading of laboratory glassware – can be equipped with 18 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Equipped with 9 A 842 nozzles and 9 A 843 nozzles
- Part of the SimpleLoad system
- Nozzle spacing 3" (78 mm) across the width, 3" (78 mm) across the depth

EAN: 4002516188971 / Material number: 11116950 / Old Material Number: 69630105D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	18
Dimensions and weight	
External dimensions, net height in inches	8 (203) in (mm)
External dimensions, net width in inches	9 3/16 (232) in (mm)
External dimensions, net depth in inches	18 15/16 (480) in (mm)
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	9 11/16 (245)
External dimensions, gross depth in inches	20 5/16 (515)
Net weight in lbs	5 (2.16)
Gross weight in lbs	7 (3.17)
Accessory pack	
A 842, injector nozzle Ø 3/16" (4 mm), length 3 9/16" (90 mm) [quantity]	9
A 843, injector nozzle Ø 3/16" (4 mm), length 7 5/16" (185 mm) [quantity]	9