


A 302/3

Module

For optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Equipped with 16 A 844 nozzles and 16 A 845 nozzles
- Part of the SimpleLoad system
- Nozzle spacing 2 3/8" across the width, 2 1/4" across the depth

EAN: 4002516188988 / Material number: 11116960 / Old Material Number: 69630203D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	32
Dimensions and weight	
External dimensions, net height in inches	5 11/16 (143) in (mm)
External dimensions, net width in inches	9 3/8 (238) in (mm)
External dimensions, net depth in inches	18 7/8 (479) in (mm)
External dimensions, gross height in inches	5 3/8 (135)
External dimensions, gross width in inches	11 7/16 (290)
External dimensions, gross depth in inches	23 1/4 (590)
Net weight in lbs	4 (1.49)
Gross weight in lbs	6 (2.67)
Accessory pack	
A 844, injector nozzle Ø 1/8" (2.5 mm), length 3 1/8" (80 mm) [quantity]	16
A 845, injector nozzle Ø 1/8" (2.5 mm), length 5" (125 mm) [quantity]	16