


A 14/1  
Lid  
for securing delicate wash items.



- Lid in stainless steel
- 1/4" x 1/4" (7 x 7 mm) perforations, 1/8" (3 mm) slat

EAN: 4002516340126 / Material number: 11131560 / Old Material Number: 69746003D

|  |                 |
|--|-----------------|
| <b>Product category</b>                        |                 |
| Large-capacity washer-disinfectors, laboratory | •               |
| Washer-disinfectors, medical                   | •               |
| Large-capacity washer-disinfectors, medical    | •               |
| <b>Product classification</b>                  |                 |
| Washer-disinfectors, laboratory                | •               |
| <b>Machine/accessory categorization</b>        |                 |
| G 7823   | •               |
| G 7824   | •               |
| G 7825   | •               |
| G 7826   | •               |
| PG 8527  | •               |
| PG 8528  | •               |
| PG 8535  | •               |
| PG 8536  | •               |
| PG 8562  | •               |
| PG 8582  | •               |
| PG 8582 CD                                     | •               |
| PG 8504  | •               |
| PG 8583  | •               |
| PG 8583 CD                                     | •               |
| PG 8592  | •               |
| PG 8593  | •               |
| PLW 8615                                       | •               |
| PLW 8616                                       | •               |
| PLW 8617                                       | •               |
| PLW 6011                                       | •               |
| PLW 6111                                       | •               |
| PLW 8505                                       | •               |
| <b>Product category</b>                        |                 |
| Cover  | •               |
| Lid  | •               |
| <b>Product characteristics</b>                 |                 |
| Material                                       | Stainless steel |
| Color  | Stainless Steel |
| <b>Application area</b>                        |                 |
| Reprocessing of laboratory glassware           | •               |
| Reprocessing of different utensils             | •               |
| <b>Dimensions and weight</b>                   |                 |
| External dimensions, net height in inches      | 18.2            |
| External dimensions, net width in inches       | 223             |
| External dimensions, net depth in inches       | 256             |
| External dimensions, gross height in inches    | 2 9/16 (65)     |
| External dimensions, gross width in inches     | 9 1/2 (240)     |
| External dimensions, gross depth in inches     | 9 7/8 (250)     |
| Net weight in lbs                              | 1 (0.20)        |
| Gross weight in lbs                            | 1 (0.35)        |