

AK 12/1  
Mesh Basket  
for holding different utensils.



- Suitable for beakers, funnels, breast milk pumps, for example
- With flexible handles
- Suitable for many upper and lower baskets as well as mobile units

EAN: 4002516343714 / Material number: 11131690 / Old Material Number: 69501204D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
PLW 8505	•
<b>Product category</b>	
Insert/module for beakers	•
Insert/module for wide-neck glassware	•
Other inserts	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	4 13/16 (122)
External dimensions, net width in inches	8 7/8 (224)
External dimensions, net depth in inches	17 1/8 (434)
External dimensions, gross height in inches	4 15/16 (125)

AK 12/1  
Mesh Basket  
for holding different utensils.



EAN: 4002516343714 / Material number: 11131690 / Old Material Number: 69501204D

External dimensions, gross width in inches	9 11/16 (245)
External dimensions, gross depth in inches	17 3/4 (450)
Net weight in lbs	2 (0.65)
Gross weight in lbs	3 (1.07)