


APLW 096

Upper basket

For the optimum loading of wide-necked laboratory glassware.



- Useable surface: 485 x 450 mm
- Vertical clearance 220 mm
- With integrated spray arm

EAN: 4002516130796 / Material number: 11149030 / Old Material Number: 69309600D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8505	•
Product category	
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Blue, Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	13 7/8 (350)
External dimensions, gross width in inches	22 1/4 (562)
External dimensions, gross depth in inches	24 1/4 (615)
Net weight in lbs	10 (4.29)
Gross weight in lbs	16 (7.09)