

TT 86

Transfer trolley

For loading and unloading washer-disinfectors and lab washers.



- Equipped with 4 lockable castors
- Simple and convenient handling thanks to an ergonomic handle
- Fitted with a removable drip tray
- Height-adjustable
- Suitable for cleaning in container and trolley washers

EAN: 4002516147725 / Material number: 11178960 / Old Material Number: 69200115D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Storage/transport system	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	911
External dimensions, net width in inches	693
External dimensions, net depth in inches	990
External dimensions, gross height in inches	44 1/2 (1130)
External dimensions, gross width in inches	31 1/2 (800)
External dimensions, gross depth in inches	47 1/4 (1200)
Net weight in lbs	49 (22.20)
Gross weight in lbs	91 (41.00)