

TT 86

Transfer trolley

For loading and unloading washer-disinfectors and lab washers.



- Equipped with 4 lockable castors
- Simple and convenient handling thanks to an ergonomic handle
- Fitted with a removable drip tray
- Height-adjustable
- Suitable for cleaning in container and trolley washers

EAN: 4002516147725 / Material number: 11178960 / Old Material Number: 69200115D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Storage/transport system	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	35 7/8 (911)
External dimensions, net width in inches	27 5/16 (693)
External dimensions, net depth in inches	39 (990)
External dimensions, gross height in inches	44 1/2 (1130)
External dimensions, gross width in inches	31 1/2 (800)
External dimensions, gross depth in inches	47 1/4 (1200)
Net weight in lbs	49 (22.20)
Gross weight in lbs	90 (41.00)