

Classic C1 Turbo Team PowerLine - SBANO
canister vacuum cleaners
with Turbobrush – ideal for pet lovers.



- High suction power – 1200 W
- Reliably removes pet hair and lint with turbo brush.
- Adjustable working height – Stainless steel telescopic tube
- Effortless vacuuming of large areas w/ 29.5 ft operating radius
- Gentle on sensitive hard floors – soft parquet brush

EAN: 4002516149064 / Material number: 11181010 /
Old Material Number: 41BAN044USA

General information	
Model number	Classic C1 Turbo Team PowerLine – SBANO
Manufacturer material number	11181010
EAN number	4002516149064
Design	
Canister vacuum cleaner with dustbag	•
Product line	Classic C1
Motor type	PowerLine
Product type	SBANO
Construction type and design	
Color	Tech blue
User convenience	
Electronic suction power setting	6-stage rotary selector
Ergonomic carrying handle	•
Ergonomic handle	•
XXL handle	•
Silence system	•
Telescopic suction wand	Comfort Compact telescopic wand
Park system for vacuuming breaks	•
Park system for storing	on one side
Mobility	
Caster wheels gentle on floors	•
Three-point maneuverability	•
Operating radius in ft (m)	29 (9)
Cord length in ft (m)	18 (5.5)
Filter system	
FilterBag type	HyClean GN
FilterBag™ volume in qt (US) (l)	4.8 (4.5) qt (US) (l)
FilterBag™ change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic FilterBag™ positioning	•
Technical data	
Power	1200
Width in inches (mm)	18 1/4 (463)
Height in inch (mm)	279
Depth in mm	8 3/4 (221)
Weight in lbs (kg)	13 (5.8)
Standard accessories	
Accessory storage	On VarioClip
Upholstery nozzle	•
Crevice nozzle	•
Dust brush	•
Parquet brush SBB 300-3	•
Turbobrush STB 305-3 (DOUBLE; LOCKED)	•

Classic C1 Turbo Team PowerLine - SBANO
canister vacuum cleaners
with Turbobrush – ideal for pet lovers.