

A 850
Holding frames
For injector module A 620.



- Arrangement: 2 x 5 fields
- For use on A 620 injector module
- For the optimum positioning of laboratory glassware
- Part of the SimpleLoad system

EAN: 4002516200703 / Material number: 11236630 / Old Material Number: 69685000D

| | |
|--|--------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Other inserts/other modules | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Color | Stainless Steel, Black |
| Number of compartments | 10 |
| Compartment width in inches | 4 1/2 (114) |
| Compartment length in inches | 4 3/8 (110) |
| Part of the SimpleLoad system | • |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Dimensions and weight | |
| External dimensions, net height in inches | 5 13/16 (147) in (mm) |
| External dimensions, net width in inches | 9 5/8 (244) in (mm) |
| External dimensions, net depth in inches | 22 13/16 (578) in (mm) |
| External dimensions, gross height in inches | 5 3/8 (135) |
| External dimensions, gross width in inches | 11 1/4 (285) |
| External dimensions, gross depth in inches | 23 1/16 (585) |
| Net weight in lbs | 2 (0.55) |
| Gross weight in lbs | 3 (0.95) |