


A 851  
Holding frames  
For injector module A 621.



- Arrangement: 3 x 7 fields
- For use on A 621 injector module
- For the optimum positioning of laboratory glassware
- Part of the SimpleLoad system

EAN: 4002516201915 / Material number: 11236640 / Old Material Number: 69685100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Other inserts/other modules	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Number of compartments	21
Compartment width in inches	2 15/16 (74)
Compartment length in inches	3 1/16 (77)
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	4 11/16 (119) in (mm)
External dimensions, net width in inches	9 5/8 (243) in (mm)
External dimensions, net depth in inches	22 7/8 (580) in (mm)
External dimensions, gross height in inches	5 3/8 (135)
External dimensions, gross width in inches	11 1/4 (285)
External dimensions, gross depth in inches	23 1/8 (586)
Net weight in lbs	2 (0.67)
Gross weight in lbs	3 (1.07)