


A 852
Holding frames
For injector module A 622.



- Arrangement: 4 x 9 fields
- For use on A 622 injector module
- For the optimum positioning of laboratory glassware
- Part of the SimpleLoad system

EAN: 4002516201922 / Material number: 11236650 / Old Material Number: 69685200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Other inserts/other modules	•
Product characteristics	
Material	Wire, rilsan-coated, Stainless steel, Plastic
Color	Stainless Steel, Grey, Black
Number of compartments	36
Compartment width in inches	2 3/8 (59)
Compartment length in inches	2 5/16 (58)
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	3 (76) in (mm)
External dimensions, net width in inches	10 1/8 (256) in (mm)
External dimensions, net depth in inches	22 3/16 (562) in (mm)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	13 (330)
External dimensions, gross depth in inches	26 13/16 (680)
Net weight in lbs	2 (0.61)
Gross weight in lbs	3 (0.99)