

A 12/2

Mesh Basket

For the gentle cleaning of delicate wash items by interrupting the spray jet.



- With perforations 1/4" x 1/4" (7 x 7 mm)
- For use in upper and lower basket
- Covers half of basket
- Not for use in PG 8536

EAN: 4002516237068 / Material number: 11238030 / Old Material Number: 69501205D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Spacer	•
Retaining frame	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	9 3/8 (237)
External dimensions, net width in inches	17 11/16 (448)
External dimensions, net depth in inches	0 9/16 (14)
External dimensions, gross height in inches	2 9/16 (65)
External dimensions, gross width in inches	10 1/16 (255)
External dimensions, gross depth in inches	18 3/4 (475)
Net weight in lbs	1 (0.39)
Gross weight in lbs	2 (0.57)