


A 11/2

Mesh Basket

For the gentle cleaning of delicate wash items by interrupting the spray jet.

EAN: 4002516340119 / Material number: 11239130 / Old Material Number: 69501103D



- With perforations 7 x 7 mm
- For use in upper and lower basket

Machine/accessory categorization	
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8505	•
Product category	
Spacer	•
Retaining frame	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	1
External dimensions, net width in inches	476
External dimensions, net depth in inches	443
External dimensions, gross height in inches	2 7/8 (70)
External dimensions, gross width in inches	20 (505)
External dimensions, gross depth in inches	17 (430)
Net weight in lbs	2 (0.71)
Gross weight in lbs	3 (1.08)