

A 320

Mesh Basket

For the gentle cleaning of delicate wash items by interrupting the spray jet.



- With perforations 1/4" x 1/4"
(7 x 7 mm)

EAN: 4002516237082 / Material number: 11295510 / Old Material Number: 69632000D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Retaining frame	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	0 7/16 (10) in (mm)
External dimensions, net width in inches	9 3/8 (237) in (mm)
External dimensions, net depth in inches	21 1/8 (536) in (mm)
External dimensions, gross height in inches	2 3/8 (60)
External dimensions, gross width in inches	9 11/16 (245)
External dimensions, gross depth in inches	22 1/16 (560)
Net weight in lbs	2 (0.47)
Gross weight in lbs	2 (0.63)