

KWT 6722 iS Built-in wine storage unit with FlexiFrame and SommelierSet for wine connoisseurs.



- Flexible storage of larger wine bottles thanks to FlexiFrame
- Simultaneous wine storage in separate temperature zones
- Odorless refrigeration section thanks to Active AirClean filters
- Present the storage of your wine in the perfect light LED lighting

EAN: 4002516235422 / Material number: 11329860 / Old Material Number: 36672200USA

•
•
•
right
•
•
LED illumination
•
XKS 3130 W
•
•
•
8
Beechwood
8
•
Yes
•
•
TouchControl
No
2
5
20
•
•
•
•
•
22 1/8 (560)
22 1/2 (570)
69 7/8 (1772)
70 1/8 (1780)
21 3/4 (550)
22 (557)
69 3/4 (1770)
21 7/8 (553)
203 (92.0)
Fixed door
Fixed door 45 (20)
Fixed door 45 (20) 8.97 (254)



KWT 6722 iS Built-in wine storage unit with FlexiFrame and SommelierSet for wine connoisseurs.




EAN: 4002516235422 / Material number: 11329860 / Old Material Number: 36672200USA

0.75 Bordeaux bottle capacity (no. of bottles)	80 bottles
Current consumption in milliampere (mA)	1000
Voltage in V	120
Fuse rating in A	10
Number of phases	1
Frequency in Hz	50-60
Wire length in ft (m)	9 (2.8)
Standard accessories	
Active AirClean filter	•



KWT 6722 iS Built-in wine storage unit with FlexiFrame and SommelierSet for wine connoisseurs.



 $\label{eq:KWT6722iS-1} \mbox{KWT6722iS-1, Interior dimensions (Installation drawings)}$