

KWT 6722 iS

Built-in wine storage unit
with FlexiFrame and SommelierSet for wine connoisseurs.



- Flexible storage of larger wine bottles thanks to FlexiFrame
- Simultaneous wine storage in separate temperature zones
- Odorless refrigeration section thanks to Active AirClean filters
- Present the storage of your wine in the perfect light – LED lighting

EAN: 4002516235422 / Material number: 11329860 /
Old Material Number: 36672200USA

Appliance category	
Wine storage unit	•
Design	
Built-in	•
Built-in integrated	•
Door hinge side	right
Convertible hinging	•
Can be installed side-by-side	•
Design	
Wine compartment lighting	LED illumination
User convenience	
Networking with Miele@home	•
Required Accessories	XKS 3130 W
DynaCool	•
SoftClose door closing damper	•
Silence system	•
Total number of shelves	8
Number of shelves	8
Shelf material	Beechwood
Number of FlexiFrame Wire Racks	8
SommelierSet	•
NoteBoard	Yes
3 D bottle presenter	•
Low-vibration storage	•
Control	
Electronic temperature display and control	TouchControl
Independent temperature control of fridge and freezer compartments	No
No. of temperature zones	2
Min. temperature setting range in °F (°C)	5
Max. temperature setting range in °F (°C)	20
Safety	
Lock function	•
Acoustic door alarm	•
Acoustic temperature alarm	•
Optical door alarm	•
Optical temperature alarm	•
Technical data	
Niche width in inches (mm), min.	22 1/8 (560)
Niche width in inches (mm), max.	22 1/2 (570)
Niche height in inches (mm), min.	69 7/8 (1772)
Niche height in inches (mm), max.	70 1/8 (1780)
Niche depth in in.(mm)	21 3/4 (550)
Appliance width in in.(mm)	22 (557)
Appliance height in in.(mm)	69 3/4 (1770)
Appliance depth in in.(mm)	21 7/8 (553)
Weight in lbs.(kg)	203 (92.0)
Type of mounting	Fixed door

KWT 6722 iS

Built-in wine storage unit

with FlexiFrame and SommelierSet for wine connoisseurs.




EAN: 4002516235422 / Material number: 11329860 /

Old Material Number: 36672200USA

Maximum fridge compartment door front weight in lbs. (kg)	45 (20)
Wine compartment in L	8.97 (254)
Total usable capacity in cu. ft. (l)	8.97 (254)
0.75 l Bordeaux bottle capacity (no. of bottles)	80 bottles
Current consumption in milliampere (mA)	1000
Voltage in V	120
Fuse rating in A	10
Number of phases	1
Frequency in Hz	50-60
Wire length in ft (m)	9 (2.8)
Standard accessories	
Active AirClean filter	•

KWT 6722 iS

Built-in wine storage unit
with FlexiFrame and SommelierSet for wine connoisseurs.



KWT6722iS, KWT6722iS-1, Interior dimensions
(Installation drawings)