

DOS K 85 flex

Dispenser module for 5-10 liter canisters

with level indicator for the automatic dispensing of liquid media.



- Dispenser monitoring – in accordance with DIN EN ISO 15883
- For dispensing liquid cleaning agent/neutralizing agent
- 13 7/8 (352) in (mm) inch suction wand for 5–10 l canisters with flexible cap
- Suction-wand conversion kit for 10–30 l canisters available

EAN: 4002516365280 / Material number: 11498910 / Old Material Number: 69747052USA

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Dispenser system	•
Product characteristics	
Production rate in ml/60 s	55
Length of suction wand in mm	13 7/8 (352)
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-240
Frequency in Hz	60
Total rated load in kW	0.01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	6 3/8 (160)
External dimensions, gross width in inches	7 7/8 (198)
External dimensions, gross depth in inches	18 7/8 (478)
Net weight in lbs	5 (2.21)
Gross weight in lbs	6 (2.43)
Standards/certificates/symbols	
CE	•