

## DOS K 85 flex

Dispenser module for 5-10 liter canisters

with level indicator for the automatic dispensing of liquid media.



- Dispenser monitoring – in accordance with DIN EN ISO 15883
- For dispensing liquid cleaning agent/neutralizing agent
- 13 7/8 (352) in (mm) inch suction wand for 5–10 l canisters with flexible cap
- Suction-wand conversion kit for 10–30 l canisters available

EAN: 4002516365280 / Material number: 11498910 /  
Old Material Number: 69747052USA

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Dispenser system	•
<b>Product characteristics</b>	
Production rate in ml/60 s	55
Length of suction wand in mm	13 7/8 (352)
<b>Standard electrical connection</b>	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-240
Frequency in Hz	60
Total rated load in kW	0.01
<b>Dimensions and weight</b>	
External dimensions, gross height in inches	6 5/16 (160)
External dimensions, gross width in inches	7 13/16 (198)
External dimensions, gross depth in inches	18 7/8 (478)
Net weight in lbs	4 (1.78)
Gross weight in lbs	5 (2.02)
<b>Standards, test certificates, and symbols</b>	
CE	•