

RF MF27/28

Adapter

for using the MF 27/28 transfer trolley with the PLW 8615/16/17 model series.



**Image
in progress**

- Retrofit release element for the MF 27/28 transfer trolley
- Installed by a service technician
- Two are required for the double-door PLW 86165 machines

EAN: 4002516352471 / Material number: 11499280 / Old Material Number: 69249201D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Adapter	•
Product characteristics	
Material	Plastic
Color	Black
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	0 (0) in (mm)
External dimensions, net width in inches	0 (0) in (mm)
External dimensions, net depth in inches	0 (0) in (mm)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	4 3/4 (120)
External dimensions, gross depth in inches	6 3/4 (170)
Net weight in lbs	1 (0.25)
Gross weight in lbs	1 (0.25)